



# interWORKS

solutions for Orthopedics

## Is your orthopedic PACS living up to your expectations?

Orthopedic surgeons face an ever-expanding set of challenges including the pressure to reduce costs and deliver faster, higher-quality patient care. To meet these challenges, Avreo designed interWORKS as a WAN optimized, web-based system, that streamlines your orthopedic imaging workflow in the office, operating room, and wherever your practice takes you. Its zero-impact installation design does not require third-party

utilities, which allows it to go where many viewers can't. Its simple, yet powerful viewer automatically displays images according to each orthopedic surgeon's preferences. Mobile reference tools provide powerful visual illustrations when discussing treatment plans with patients. In addition, interWORKS provides an unlimited number of licenses allowing your facility to grow with predictable monthly costs.

## Orthopedic Workflow Solutions



## What Do I Need to Know?

**Our organization uses a reading group to provide interpretation for specific studies. Can you send these studies to their PACS?**

Yes. A technologist can send these studies as needed to your reading group. In addition, your reading group can read remotely from interWORKS if needed.

**Can interWORKS receive orders from my Practice Management/EHR?**

Yes. To maximize workflow efficiency, interWORKS can receive orders from your practice management/EHR system via HL7, which minimizes duplicate data entry requirements. Our visual DayPlanner then displays these orders to help technologists manage their schedule within a single application.

**Do you support orthopedic templating tools?**

Yes. Avreo integrates with orthopedic templating software to enable surgeons to plan and template orthopedic procedures and reduce or eliminate their need for printed film. Templated images and the planning report are available in interWORKS for future reference.

## What Does interWORKS Include?

- Cross-platform support including Windows®, Mac OS® X via Parallels, Apple® iOS, and Android™
- Individualized workflow for each orthopedic surgeon in your practice.
- Case notes and scanned documents to manage the unique workflow needs of orthopedic practices.

# interWORKS solutions for

## Orthopedics

### Recommended Configuration

- interWORKS VM hosted in your organization's data center.


### Why interWORKS?

- Easily accepts imaging orders from your PM/EHR, and provides DICOM modality worklists minimizing data entry and errors. Close the loop by providing hyperlink access to orthopedic imaging from your EHR.
- Empower orthopedic surgeons to read virtually anywhere with WAN optimized, web-based image access.
- Status driven workflow configurable to the needs of orthopedic practices.
- Quickly document case notes via digital dictation or quick-pick templates.
- Storage agnostic, vendor neutral archiving.
- Standards based integration via DICOM and HL7.

### Orthopedics Highlights

- Robust tools including distance, angle and Cobb angle measurements, spine labeling, 3D localization, window/level, zoom, and magnify.
- Configurable workflow and preferences per orthopedic surgeon. This includes layout preferences, worklists, and rules based auto-loading of comparison studies.
- Import and export DICOM CDs and include embedded viewer for ease of access.
- Coordinate care with referring physicians by providing direct access to the application or mobile app.

## Easily Access Orthopedic Imaging throughout your Enterprise



## Let's Start a Conversation!

We want to discuss the current challenges that you are facing so we can fully understand the needs of your organization and offer solutions to alleviate these issues.

### To Learn More

 [www.avreo.com](http://www.avreo.com)

 1-866-286-8082

 Demonstrations

